


AMPOL

From Poland, for Europe.

XporT

NMR Sample Tubes
Probówki-NMR
2010




CHOOSE THE BEST...

... CHOOSE AMPOL



AMPOL PRODUCES ITS NMR TUBES EXCLUSIVELY IN THE USA. OUR QUALITY STANDARDS ARE VERY HIGH AND OUR PRICES ARE COMPETITIVE. YOUR SATISFACTION IS UNCONDITIONALLY MONEY-BACK GUARANTEED!



AMPOL, P.O. Box Nr. 60 18-404 Lomza, Poland
Tel. 48-86-218-9908; Fax. 48-86-212-5034

For your NMR Equipment

For more details,

visit us at www.ampolnmr.com

HIGH CONCENTRICITY NMR SAMPLE TUBES

from AMPOL... YOUR BEST VALUE!

Ampol maintains strict concentricity and very high tolerances in its manufacture of Precision NMR Sample Tubes. Once you begin using Ampol NMR Tubes, you may discover that shimming is not required when you change the tubes in your

NMR probe. Please note that our most popular

AP5-400-7 NMR tube, recommended for 300/400MHz NMR equipment comes with 15 micron tolerance for both, concentricity and camber. Our AP5-600-7, second most popular NMR tube, is recommended for 500/700 MHz NMR units and retains tolerance of 10 microns.

AMPOL'S 2010 NMR SAMPLE TUBE PRICE LIST

AMPOL, P.O. Box Nr. 60 18-404 Lomza, Poland Tel. 48-86-218-9908; Fax. 48-86-212-5034

To establish permanent order record, please expedite your orders by fax, e-mail or by post.

DISCOUNTS are granted on purchases in bulk. Save € 3,500 on purchase of € 10,000 or more. This rate also applies to assorted purchases.

Discount rate:	10%	15%	20%	25%	30%	35%
Total purchase required:	500	1,000	2,500	5,000	7,500	10,000

Order 10,000 AP5-60-7 (high volume use NMR tube) and you pay only € 0,8694/tube

5 mm NMR Sample Tubes. Tube end caps, No. TEC-5, supplied in lots of 100(€ 6,42) and 1000(€ 48,15).

Colors available: red, blue, green, yellow, black. All 5-pack-tubes come with red caps.

NMR Sample Tube No.	MHz NMR Rating	Tube Length cm	Tube Length inch	Wall thickness mm	Concentricity microns	Camber microns	Price per Tube EURO	Package Size (Tubes/Pack)
AP5-1000-7	900	17,8	7	0,385	2	2	33,17	5 tubes
AP5-1000-8	900	20,3	8	0,385	2	2	37,45	5 tubes
AP5-900-7	900	17,8	7	0,385	3	3	19,26	5 tubes
AP5-900-8	900	20,3	8	0,385	3	3	22,47	5 tubes
AP5-800-7	800	17,8	7	0,385	5	5	14,98	5 tubes
AP5-800-8	800	20,3	8	0,385	5	5	17,12	5 tubes
AP5-600-7	500/700	17,8	7	0,385	10	10	8,56	5 tubes
AP5-600-8	500/700	20,3	8	0,385	10	10	9,63	5 tubes
AP5-400-7	300/400	17,8	7	0,385	15	15	6,42	25 tubes
AP5-400-8	300/400	20,3	8	0,385	15	15	7,49	25 tubes
AP5-200-7	100/300	17,8	7	0,385	20	20	4,28	25 tubes
AP5-200-8	100/300	20,3	8	0,385	20	20	4,92	25 tubes
AP5-100-7	60/100	17,8	7	0,385	30	30	2,03	100 tubes
AP5-100-8	60/100	20,3	8	0,385	30	30	2,25	100 tubes
AP5- 60-7	60	17,8	7	0,385	40	40	1,28	100 tubes
AP5- 60-8	60	20,3	8	0,385	40	40	1,50	100 tubes

10mm NMR Sample Tubes, supplied with end caps No. TEC-10

AP10-600-7	400/600	17,8	7	0,60	10	10	32,10	5 tubes
AP10-600-8	400/600	20,3	8	0,60	10	10	40,66	5 tubes
AP10-300-7	200/300	17,8	7	0,60	20	20	26,75	5 tubes
AP10-300-8	200/300	20,3	8	0,60	20	20	31,03	5 tubes
AP10-200-7	60/200	17,8	7	0,60	40	40	14,98	5 tubes
AP10-200-8	60/200	20,3	8	0,60	40	40	17,12	5 tubes

AMPOL
Box Nr. 60
18-404 LOMZA-6, POLAND
TEL. 48-86-218-9908
FAX. 48-86-212-5034
www.ampolnmr.com

